



Hill-Rom VersaCare Bed - Fully Refurbished

The VersaCare bed is a full-featured patient care bed system that offers the lowest height in acute care, an integrated air surface that helps automate the delivery of wound care therapy, and many features to reduce the risk of caregiver injury. In addition, this latest innovation from Hill-Rom offers superior patient comfort and intuitive controls to help caregivers effectively and easily operate the system to its fullest capacity.

Features:

- Built-in scare and three-level bed exist system
- Low height and low chair position
- Docking patient pendant controls
- Dedicated heel relief zone
- SafeView alerts
- Cord wrap clips and IV pole storage
- HandsFree foot controls
- Nightlight
- FlexAFoot bed length adjustment
- Lowest height of any acute care bed helps reduce the risk and impact of patient falls.
- Firm, supported perimeter of surface enables safer patient ingress and egress.
- Minimized siderail gaps reduce opportunities for entrapment. • Help patients achieve greater comfort in LowChair® position
- Surface thicknesses ranging from seven to eight inches surpass other surface thicknesses for patient comfort.
- Bed height: low position - 18" and high position - 37"
- Overall length: extended - 94.5" and retracted - 82.5"
- Width: side rails stowed - 37" and side rails up - 40"



Why an Auxo Medical Refurbished Bed or Stretcher?

Auxo beds and stretchers are refurbished to the some of the highest standards in the industry. We completely disassemble each bed and stretcher, and completely rework or replace all worn parts and return each bed and stretcher to OEM standards. Each bed and stretcher is reassembled, tested and is warranted against defects.