



# auxoMEDICAL

## Ritter 75E Procedure Chair - Fully Refurbished

The Ritter 75 Evolution surgical chair has the flexibility you need for a variety of uses including: Family and General Practice, OB/GYN, Internal Medicine, Dermatology, Oral, Plastic & General Surgery.

### Features:

- Minimum Height - 26" (66cm)
- Maximum Height - 42" (106.7cm)
- Width (top) - 27" (68.6cm)
- Weight - 475 lbs.
- Tilt - Trendelenburg to 45 Degrees
- Foot Control - To power height, tilt, back & foot
- Electrical - 120 VAC, 60Hz, 12AMP, UL
- Standard - Power adjustment, drain pan, articulating head rest, adjustable stirrups, accessory rails, foot section lockable in four procto positions, paper roll holder, paper straps, quick release upholstery tops.

### Accessories:

- Knee crutches (standard or articulating)
- Chair Arms
- Castor base (standard or articulating)
- Facial pad (fits standard and special procedures headrest)
- Restraint straps
- Base rails (for seat section)
- IV pole with clamp
- Armboard with pad and clamp (fixed or articulating)
- Hand surgery armboard with pad, leg and clamp
- Special procedures headrest
- Wrist support headrest
- Welch-Allyn light bracket kit
- Side rails



### Why an Auxo Medical Refurbished Table or Chair?

Auxo Medical tables are refurbished to the highest standards in the industry. We disassemble each table and chair, make necessary cosmetic upgrades, replace hydraulic components and seals and return the table to original OEM Specifications. We do not substitute cheaper non-OEM parts like many of our competitors. Each table is reassembled, tested and is warranted against defects.