



auxoMEDICAL

Steris 3085 OR Table - Fully Refurbished

The AMSCO 3085 General Surgical Table is one of the most reliable and popular surgical tables. Trusted for years by healthcare facilities across the world, the 3085 is a cost-effective solution for the general OR.

Features:

- Return-to-level button enhances patient safety
- Reverse orient button simplifies patient positioning
- Radiolucent Orthopedic Extension quickly converts this workhorse general surgical table into a fracture table
- Bariatric Width Extenders stretch the table width to 28" to accommodate the girth of up to a 1,000 lb patient
- Over 300 surgical table accessories are available to accommodate virtually any surgical specialty

Specifications:

- Table Top Length: 76"
- Table Top Width: 20" (28" with extenders)
- Table Height Range: 27" to 44"
- Patient Weight Capacity: 1,000 lbs. (454 kg) normal orientation
- Hand control features: normal/reverse orientation button + return to level button
- Trendelenburg / Reverse: 25° / 25°
- Lateral Tilt (left/right): 18° / 18°
- Head Section: +90° to -90°
- Back Section: +55° to -25° (normal) and +80° to -105° (reverse)
- Leg Section: +80° to -105° (normal) and +55° to -25° (reverse)
- Flex/Reflex: 20° / 25°



Accessories Available: Arm Boards & Hand Tables, Foot Control, X-Ray Boards, Patient Straps & Restraints

Why an Auxo Medical refurbished Surgical Table?

Auxo Medical tables are refurbished to the highest standards in the industry. We disassemble each table, make necessary cosmetic upgrades, replace hydraulic components and seals and return the table to original OEM Specifications. We do not substitute cheaper non-OEM parts like many of our competitors. Each table is reassembled, tested and is warranted against defects.