



Zimmer ATS 2000 Tourniquet System - Refurbished

The Zimmer ATS 2000 Automatic pneumatic Tourniquet System's built-in safety, reliability, and convenience frees you to focus on surgery. Dual displays that monitor actual cuff/bladder pressure, precise microprocessor controls, and 1-hour battery backup help assure patient safety. User-set alarms warn of unacceptable or hazardous conditions.

Features:

- Dual controls, dual hosing, dual cuffs, and dual readouts provide independent and precise monitoring and performance for bilateral or Bier Block procedures
- Positive-locking connectors supply safer, leak-free inflation
- Angled ports keep connectors clear of the surgical site
- Visual and adjustable audio alarms readily detect pressure variations
- Easy-to-read display adjusts easily
- Simple setup includes user-defined default pressure and time settings
- Lightweight portable unit can be mounted on an I.V. pole or placed on a table
- One-hour battery backup – just in case

Specifications:

- 100-240V - 50/60 Hz - 120 VA
- Adjustable Pressure: 1mmhg - 45mmhg
- Adjustable Time: 1 minute - 240 minutes
- Power Cord
- Height: 12.75 in.
- Width: 10 in.
- Depth: 10.5 in.
- Weight: 19.8 lbs
- * Cuffs Not Included

Accessories: Cuffs, IV Pole with Basket



Why purchase Auxo Medical Reconditioned Biomedical Equipment?

Auxo thoroughly inspects and re-certifies all of its reconditioned biomedical equipment. Batteries are checked or replaced, and each unit is cosmetically repaired. Calibrations are made and all safety inspections will be completed. In most cases our equipment warranties meets or exceeds those of the OEM.