

EDAN



M3B Vital Signs Monitor

With SpO₂ (Nellcor Oximax™ optional) and Respironics CO₂ monitoring technologies, the M3B ensures effective capnography monitoring for mechanically ventilated and non-intubated patients during continuous long-term monitoring.

Measured Parameters: SpO₂ + CO₂



72 h
Trend Review

800
Alarm Review



Care for Health

M3B

Vital Signs Monitor

- 5.6" high resolution color TFT-LCD screen
- Lightweight & portable design
- User-friendly interface
- Flexible configurations to meet different clinical needs
- Real-time measurements display
- Trend table display
- Built-in rechargeable Lithium-ion battery for 8 hours of continuous work
- USB flash disk support
- Barcode scanner support
- Bi-directional with MFM-CMS central monitoring system
- Compatible with PatientCare Viewer PC data management software
- Wired LAN support
- HL7 support via XML files or MFM-CMS
- Nurse call
- Suitable for adult, pediatric and neonatal patients

Seamless Connectivity




Mounting Solutions



Wall Mount



Rolling Stand



Bed Rail Clamp

Specification

Configuration	SpO ₂	Respiration CO ₂ (Sidestream/Mainstream)
M3B(SpO ₂ + CO ₂)	✓	✓
M3B(CO ₂ only)	—	✓

EDAN Edan Instruments, Inc.

■ 3/F – B, Nanshan Medical Equipment Park, Nanshan Rd 1019#, Shekou, Nanshan Shenzhen, 518067 P.R. China
 ■ Tel: +86-755-26898326 ■ Fax: +86-755-26898330 ■ www.edan.com.cn ■ Email: info@edan.com.cn

■ All rights reserved. Features and specifications are subject to change without notice



CE 0123

US-PM-M3B-V1.2-20130514