



auxoMEDICAL

Steris Harmony LA300 - Fully Refurbished

Lighting designed for patient and caregiver with low profile, high quality illumination for critical care and labor and delivery settings. The Harmony LA300 Examination Lighting System is designed specifically for minor procedures in an ICU or ER. The LA300 provides high-intensity, color-corrected lighting and its compact mobile, ceiling- or wall-mount options can be installed anywhere procedure lighting is necessary.

Features:

- Small and extremely maneuverable
- Cool, color-corrected light
- Future-ready design protects your investment.
- Modular, flexible design easily adapts as your needs change.
- All lightheads are camera-ready.
- Intensity level is adjustable from the bezel
- Intensity: 24,000 to 60,000 lux
- CRI of 94
- Color temperature of 4400K
- Configurations include wall-mount, ceiling-mount, or mobile stand
- Focal range: 34 - 59"
- Pattern size: 9" (229 mm)
- Mobile base size: 22¼" x 22" (565 mm x 558 mm)
- Mobile base height, weight: 69" (1732 mm), 76 lb.(35 kg)
- Horizontal extension arm: 31½" (800 mm)
- Spring arm: 40" (1016 mm)



Why purchase Auxo Medical Refurbished Lights?

Auxo Lights, halogen and LED, are refurbished to the highest standards in the industry. We completely disassemble each light and arm assembly, and completely rebuild or replace all worn parts and return each assembly to OEM standards. Each Light is reassembled, tested and is warranted against defects.