

## Steris Harmony LA500 Lights - Fully Refurbished

---

The Steris Amsco Harmony LA500 surgical lighting system is a great choice for almost any operating room. The Harmony delivers a high illumination intensity of up to 14,780 foot candles (160,000 Lux) by using the 7 step intensity control, either on the light head or the wall mounted intensity controller. The intensity controller is part of a control center that has a modern LCD display. The backup lighting system means the surgeon will still have the light needed to continue, even if the main bulb burns out.

### Features:

- Control Center with Light Intensity Control and Diagnostics
- Available configurations include single, dual, or triple heads or with a combination of light heads with cameras, monitor arms, etc.
- Wavelens focusing technology
- Light Intensity Controls on Heads and Control Center
- Heat Filtering with Cold Mirror Reflectors
- Shadow Reducing Design
- Color Corrected
- Lightweight Suspension System
- Automatic Lamp Change Over System
- Standard set is configured with Dual Heads
- Intensity: 14,870 Foot Candle (160,000 Lux)
- Pattern Size: 6" - 11"
- Depth of Field: 42" (small pattern)
- Focal Range: 32" - 53" (small pattern)
- Color Temp: 4400K
- CRI: 94
- Electrical: 120V / 5 Amps



### Why purchase Auxo Medical Refurbished Lights?

Auxo Lights, halogen and LED, are refurbished to the highest standards in the industry. We completely disassemble each light and arm assembly, and completely rebuild or replace all worn parts and return each assembly to OEM standards. Each Light is reassembled, tested and is warranted against defects.