



Amsco 2080 Manual OR Table - Fully Refurbished

The Amsco 2080M manual surgical table is a multi-purpose operating table with the reliability and simplicity of manual operation. This workhorse is still very popular for its functionality and value and continues to provide a superior solution for specific clinical needs. The 2080 provides a full range of articulation, meeting the needs of virtually every general surgical procedure. The addition of the x-ray top mounted on the table top facilitates taking x-rays during the procedure. Despite being a manual operated table the controls are very easy to use and extremely intuitive.

Features:

- Pedestal includes tabletop lift cylinder and support column with bearing-mounted saddle frame fully enclosed by stainless-steel, telescoping shrouds.
- Each shroud is one-piece construction
- 300 lb weight capacity,
- Side rails for accessories
- Adjustable head rest.
- Back up/down
- Leg up/down
- Trendelenberg/Reverse Trendelenberg
- 20-inch wide tabletop divided into four hinged sections
- Cast aluminum sections



Specs:

Table Height: 45"-27"

Trendelenberg and Reverse Trendelenberg: 20°

Lateral Tilt: ±10°

Length: 72 1/4"

Width: 20"

Width w/Rail: 21 5/8"

Why an Auxo Medical Refurbished Table or Stretcher?

Auxo tables and stretchers are refurbished to the some of the highest standards in the industry.

We completely disassemble each table and stretcher, and completely rework or replace all worn parts and return each table and stretcher to OEM standards. Each table and stretcher is reassembled, tested and is warranted against defects.