



# auxoMEDICAL

## Abbott Plum A+- Fully Refurbished

The Abbott® Plum A+® Infusion Pump is a general-purpose infusion device that provides precise delivery of multiple therapies across the general spectrum of clinical care. The Plum A+® has 3 separate channels with 3 individual lines to the patient. Features fully automated secondary drug delivery, concurrent flow, preferred dosing and delivery rate (15 options), variable occlusion pressure, and bar code reader.

### Features:

- Deliver infusions with confidence with the unique PlumSet™ technology
- Automated second-line delivery — eliminates the need to raise and lower infusion containers, improving ease-of-use and efficiency
- Concurrent delivery of 2 medications through a single patient line
- Air-in-Line detection and elimination — air trap captures up to 2 mLs of air
- Accumulated air can be removed by automated back-priming without disconnecting the patient line, reducing the risk of contamination
- Deliver secondary infusions via IV container or syringe
- Reliable performance, steady, consistent and proven
- One integrated system works across the entire spectrum of clinical care
- Single channel (2 concurrent medications) or triple channel (6 concurrent medications)
- Choose from a variety of programming options:
  - Automated piggyback and concurrent delivery
  - Programmable standby settings
  - Multistep delivery
  - Loading dose automation
  - Programmable delayed starts
- Dimensions: Approx. 8"(h) x 8"(w) x 6"(d)
- Weight: Approx. 9.5 lbs



### Why purchase Auxo Medical Reconditioned Biomedical Equipment?

Auxo thoroughly inspects and re-certifies all of its reconditioned biomedical equipment. Batteries are checked or replaced, and each unit is cosmetically repaired. Calibrations are made and all safety inspections will be completed. In most cases our equipment warranties meets or exceeds those of the OEM.