



Biolight A5 AnyView Modular - Fully Refurbished

The Biolight A5 patient monitor is a versatile patient monitor, capable of being utilized in virtually any area of your healthcare facility. The large 13.75" TFT Color Display makes viewing patient data easy. The Biolight A5 comes standard with the EMS1 parameter module, which includes the following parameters: ECG (3 or 5 lead), BLT Digital SpO2, NIBP, Pulse Rate, Temp, Invasive BP (x2)

Features:

- Global-Leading all-in-one module with display
- Built-in battery and thermal printer
- Standard model includes EMS1 Module with BLT Digital SpO2, Pulse Rate, ECG (3 or 5 Lead), NIBP
- Large, Anti-glare color TFT-LCD display
- Powerful measurement, ergonomic and flexible design
- Fanless cooling system, which keeps the monitor working quietly
- Dozens of parameters available with a variety of modules assembled as needed
- Easily upgradeable, saving potential monitor replacement cost
- Advanced modular monitoring
- Connects to independent large display for multi-display
- Long life Lithium battery
- 4 Module Slots
- 4 USB ports (peripheral devices interface)



Specs:

- Size: 350mm x 340mm x 207mm
- Weight: <6.5kg
- Power Voltage: AC 100-240V 50/60Hz
- Power Input: <150VA
- Input Current: 1.7~0.8A
- 12.1" Color Anti-glare TFT-LCD
- Resolution: 800x600 pixels

Why purchase Auxo Medical Reconditioned Biomedical Equipment?

Auxo thoroughly inspects and re-certifies all of its reconditioned biomedical equipment. Batteries are checked or replaced, and each unit is cosmetically repaired. Calibrations are made and all safety inspections will be completed. In most cases our equipment warranties meet or exceed those of the OEM.