



## BladderScan BVI 9400 - Fully Refurbished

The BladderScan BVI 9400 bladder volume instrument, with NeuralHarmonics technology, is a portable 3D ultrasound device that quickly, accurately, and non-invasively measures urinary bladder volume and post-void residual (PVR) in adults and small children.

### Features:

- Quick and accurate measurement of urinary bladder volume anytime, especially post-void residual (PVR)
- Distinct scan modes for men, women, and small children
- Precision aiming via console and probe, aiming assistance helps staff to learn and use
- Customizable savings calculator tracks instrument usage and UTI costs saved
- Onboard printer for patient records or reimbursement
- Measures bladder volume accurately and non-invasively on a wide range of patients
- Helps assess urinary retention and postoperative urinary retention (POUR)
- Helps evaluate many common urological problems (LUTS)
- Helps prevent unnecessary catheterization
- Helps reduce rates of catheter-associated urinary tract infection (CAUTI)
- Helps differentiate between types of incontinence to determine appropriate care
- Helps improve efficiency of healthcare professionals by reducing costs and saving staff time



### Why purchase Auxo Medical Reconditioned Biomedical Equipment?

Auxo thoroughly inspects and re-certifies all of its reconditioned biomedical equipment. Batteries are checked or replaced, and each unit is cosmetically repaired. Calibrations are made and all safety inspections will be completed. In most cases our equipment warranties meets or exceeds those of the OEM.