



Clinton Bariatric Power Table - Fully Refurbished

The Clinton Open Base, Extra Wide Bariatric Power Table with Adjustable Backrest offers ultra strong construction and a user-friendly mechanism.

Features:

- Extra wide, bariatric model with choice of two widths
- Dual-cylinder, pneumatic backrest is adjustable from either side and raises 80 degrees from flat
- Easy-clean ABS plastic enclosed lift
- Electric power height adjustment with foot control
- Ultra-strong, uni-frame base and top construction
- Heavy duty, all welded, lift mechanism
- Durable gray powder-coated steel frame
- Self-lubricating pivot points
- Smooth 110 volt electric height adjustment
- Power supply and motor meets UL 60601-1 and CAN/CSA C22.2 601.1
- Hands-free low voltage foot control
- Easy service, self contained, electrical components
- Seamless, rounded corner, upholstered top
- 2" (5 cm) ultra-firm padding
- Low height may help achieve ADA tax credit
- Customize with options



Specs:

Dimensions: 72"(L) x 34"(W) x 18" – 35"(H)

Backrest Dimensions: 38" (L) x 34"(W)

Weight Limit: 550 lbs. (249.47 kg)

Why an Auxo Medical Refurbished Table or Stretcher?

Auxo tables and stretchers are refurbished to the some of the highest standards in the industry.

We completely disassemble each table and stretcher, and completely rework or replace all worn parts and return each table and stretcher to OEM standards. Each table and stretcher is reassembled, tested and is warranted against defects.