



ConMed Excalibur Plus PC- Fully Refurbished

The ConMed™ Excalibur Plus PC™ is the most powerful isolated ESU available in its class. It allows the surgeon to adjust the power settings from the surgical field with a standard ES pencil. The Excalibur Plus PC™ patented Controlled Output Voltage in CUT maintains minimum necessary voltage for smooth cutting through various impedance tissues with minimal thermal damage. The Aspen Return Monitor (ARMU) is the most sensitive pad monitoring system available. It combines automatic feedback circuitry and visual displays to provide the most sensitive contact quality monitoring system available.

Features:

- Most powerful isolated ESU available.
- Allows the Surgeon to adjust power settings from the surgical field with a standard ES pencil.
- Excalibur Plus PC Patented Controlled Output Voltage in CUT maintains minimum necessary voltage for smooth cutting through various impedance tissues with minimal thermal damage.
- Electrosurgical Generator with Cut, Coag, Bipolar and Blend.
- Includes a Mono-Polar Foot Pedal (Bi-Polar Foot Pedal is Optional Below).
- 300 watts of output
- The Excalibur Plus PC offers a scalpel like cut and three blend modes.
- Pure Cut Maximum Peak to Peak Voltage is 2000 Volts
- Blend 1- Maximum Peak to Peak Voltage is 2000 Volts
- Blend 2- Peak Voltage of 2200 Volts
- Blend 3- The Excalibur plus PC comes equipped with the maximum Peak to Peak Voltage of 2400 Volts.
- Weight: 40 lbs (18.5 kg)
- Dimensions: 6.25" H x 18" W x 22.75" L (15.9 x 45.8x 57.8 cm)



Why purchase Auxo Medical Reconditioned Biomedical Equipment?

Auxo thoroughly inspects and re-certifies all of its reconditioned biomedical equipment. Batteries are checked or replaced, and each unit is cosmetically repaired. Calibrations are made and all safety inspections will be completed. In most cases our equipment warranties meets or exceeds those of the OEM.