



auxoMEDICAL

ConMed Sabre 2400- Fully Refurbished

The ConMed™ Sabre 2400™ is a full-power electrosurgical generator to accommodate the needs of all surgical procedures. The Sabre 2400™ easily meets the high demanding performance needs of electrosurgery. Features independent Monopolar outputs & enhanced bipolar performance. It is safe, convenient, cost-effective, and can be programmed for standardization.

Features:

- High Power Electrosurgical Unit with Mono-Polar Foot Pedal (Bi-Polar Foot Pedal is Optional Below).
- Full-Power Electrosurgical Generator to accommodate the needs of all Surgical Procedures.
- Features Independent Monopolar Outputs, Enhanced Bipolar Performance.
- Is Able to be Programmed for Standardization.
- Combines Automatic Feedback Circuitry and Visual Displays to provide the most Sensitive Contact Quality Monitoring System Available Today.
- Continuously Monitors the Safety of the Monitoring Patient to Ground System Pad Interface with a Dual Foil Ground Pad.
- Alerts the Staff if a Pad Fault Condition occurs and Disarms the ESU.
- ARM Operates in a Preset Safety Range of 10-150 ohms.
- The Assist Nar Graph provides a Visual Alarm Set Point and Eliminates the Possibility of Pad Site Burns Attributable to Machine Error.



Specs:

- Weight: 32 lbs (14.5 kg)
- Dimensions: 6.75" H x 14.25" W x 21.5" L (17.1 x 36.2 x 54.6 cm)
- Power Decrease: 1% Typical. Pure Cut Set to 150 watts for 20 seconds
- Input Voltage: 104-127 V (rms), 50/60 Hz, Calibrated at 120 V
- Input Current: Pure Cut at 300: 7.0A Max

Why purchase Auxo Medical Reconditioned Biomedical Equipment?

Auxo thoroughly inspects and re-certifies all of its reconditioned biomedical equipment. Batteries are checked or replaced, and each unit is cosmetically repaired. Calibrations are made and all safety inspections will be completed. In most cases our equipment warranties meets or exceeds those of the OEM.