



## Cox 6000 Rapid Dry Heat - Fully Refurbished

The COX Rapid Heat Transfer Sterilizer is a fast, simple solution to your sterilization needs. Because it utilizes heated air that moves at a rate of 2,500 feet/minute, it is able to sterilize objects in as little as six (6) minutes — faster and easier than any other sterilizer available. No longer will you and your staff be bothered with all the problems inherent with other forms of sterilizations. Since the dry heat system won't corrode or rust your instruments, you won't have to replace them so often — thus saving you money.

### Features:

- Promotes longer instrument life (Does not corrode, pit, or crack)
- Eliminates need for water, rust inhibitors or chemicals
- Faster instrument turnaround – Reduces instrument investment
- Lower maintenance cost – No toxic residue or effluent

### Sterilization Cycles:

- Cycle 1: 6 minutes for unwrapped instruments
- Cycle 2: 8 minutes for hand pieces and medical drills
- Cycle 3: 12 minutes for wrapped instruments

### Included with Sterilizer:

- One Drawer
- One Removable Instrument Basket
- Basket Removal Tool
- Cooling Rack
- Instrument Rack (for wrapped instruments)



### Why purchase Auxo Medical Reconditioned Biomedical Equipment?

Auxo thoroughly inspects and re-certifies all of its reconditioned biomedical equipment. Batteries are checked or replaced, and each unit is cosmetically repaired. Calibrations are made and all safety inspections will be completed. In most cases our equipment warranties meets or exceeds those of the OEM.