



DRE Lucerne - Fully Refurbished

With the DRE Lucerne 360, you get premium features for as much as 30% less than the cost of competing surgery tables. Crafted for quality, the Lucerne 360 provides intuitive movement and stability, setting new standards in patient positioning. It features a C-Arm compatible table top with a 225 kg (500 lb) lift capacity. The Lucerne 360 has many advanced positioning options, such as kidney bridge positioning, 360° horizontal rotation (180° each direction) and a horizontal slide motion that helps create a greater field of view for C-arm use.

Features:

- C-arm interface capabilities to meet imaging requirements
- Remote control provides high operation safety and precise patient positioning with a safety lock
- Silent electro-hydraulic movements
- Mobile anti-static castors with a floor lock
- Easy access to electro-hydraulic system for service
- 500 lb. (225 kg) lift capacity
- Kidney bridge and flex positioning
- Positioning capabilities include:
 - Full 360° horizontal rotation
 - Rotation 0 ~ 180° each direction
 - Reverse Trendelenburg anteverision 30°/ retroversion³30°
 - Back Plate Adjustment +75° ~ -20°
 - Head Plate Adjustment 40°/ 90°
 - Leg Plate Adjustment: 90°
 - 4 Leg Plate Angle 0 ~ 90°
 - 105° leg decline
 - Longitudinal Shift



Why an Auxo Medical Refurbished Table or Stretcher?

Auxo tables and stretchers are refurbished to the some of the highest standards in the industry.

We completely disassemble each table and stretcher, and completely rework or replace all worn parts and return each table and stretcher to OEM standards. Each table and stretcher is reassembled, tested and is warranted against defects.