




## Datascope Passport 2 - Fully Refurbished

The Datascope Passport 2 Patient Monitor makes monitoring simple. It features an incredible intelligent display interface which allows the user to evaluate patient status at a glance. Data is presented clearly, concisely and logically. Waveforms are aligned to digital data for immediate interpretation. Graph and list displays are easily configured to personal preference. The 4 Auto-Set Alarms program is quickly configured and extremely easy to use. Three waveforms will automatically reconfigure to six when needed. Colorful numerics automatically re-adjust to maximum size for optimal display visibility.

### Features:

- Dedicated function keys provide one-touch access to all vital functions
- Colorful waveforms with automatic display configuration
- Six waveforms
- Auto-adjustable large numerics
- 3 or 5-lead ECG
- Lead-selectable respiration
- Non-invasive blood pressure,
- Temperature
- Masimo SET motion tolerant SpO2
- Built-in PCMCIA ports allow you to upgrade software, enhance memory, and transfer data
- Integrated power supply with lithium-ion battery technology
- BP cuff and hose, 1 Spo2 Finger sensor and cable,
- ECG cable and leads and recorder paper.
- Weight: 13 Lbs
- Auto-set alarm function
- 6 minute OXY-CRG trend for neonatal and pediatric applications graph and list trending
- Two trace thermal Printer



\*Mounting equipment sold separately.

### Why purchase Auxo Medical Reconditioned Biomedical Equipment?

Auxo thoroughly inspects and re-certifies all of its reconditioned biomedical equipment. Batteries are checked or replaced, and each unit is cosmetically repaired. Calibrations are made and all safety inspections will be completed. In most cases our equipment warranties meets or exceeds those of the OEM.