



Draeger Fabius GS Anesthesia Machine - Fully Refurbished

In an environment increasingly focused on cost issues, the anesthesia workplace Fabius GS offers an unprecedented price-performance ratio that will reduce your total cost of ownership. The ergonomic design of the Fabius GS facilitates efficient use and contributes to a positive, productive anesthesia workplace. Additionally, its modular concept allows you to configure precisely the workstation you require.

Features:

- High performance ventilation
- An integrated, high contrast color monitor
- Highly maneuverable
- Convenient central brake
- open, modular architecture
- Low-flow ventilation
- Standard mounting rails
- Medical-grade power outlet strips
- Easy Integration into your existing hospital information system
- 2 vaporizer system
- Ventilation modes: Pressure control, volume control, Pressure support (upgrade), w/ electronic PEEP
- 3 gas system: O₂, N₂O, Air / O₂ & N₂O Tank Yokes
- Aux O₂ Flowmeter
- AGSS Scavenging system
- CLIC Soda Sorb System
- Weight (base unit without vaporizers or cylinders): 296 lbs (134.2 kg)
- Dimensions: 30" (W) × 52"(H) × 33" (D)
- Power Supply: 100-240 VAC, 50/60 Hz, 2.3 A max.
- Battery (supports ventilator and monitor): > 45 min



Why purchase Auxo Medical Reconditioned Biomedical Equipment?

Auxo thoroughly inspects and re-certifies all of its reconditioned biomedical equipment. Batteries are checked or replaced, and each unit is cosmetically repaired. Calibrations are made and all safety inspections will be completed. In most cases our equipment warranties meets or exceeds those of the OEM.