



Dual Bay Economy Sink- IR Operated Water - Fully Refurbished

The ES-Series Surgical Scrub Sinks were developed to answer space and economical constraints where needed. These sinks have smaller outside dimensions but still maintain a deep basin to allow for proper scrubbing. This surgical scrub sink includes; infrared sensor operated water, under counter mixing valve, swivel gooseneck spout with aerator, deep sloping basin, flow control valve, 24 volt plug-in transformer. Unit will be supported by a Z-bracket and a main support bracket supplied with the sink.

Features:

- ES47 has 47" wide sink
- Deep sloping basin
- Adjustable thermostatic controlled pressure regulating mixing valve, one per station mounted under sink
- High mount swivel gooseneck spout with adjustable swivel spray aerator
- Flow control valve
- Backflow preventers
- Infrared sinks are equipped with a 24 volt plug in transformer
- Z-brackets (See specifications on rear panel for in-wall support)
- 14 Gauge 304 stainless steel polished and ground to a #4 finish
- Dimensions: 47"W x 20"D x 43"H
- Net Weight 110 Lb



Why purchase Auxo Medical Reconditioned Biomedical Equipment?

Auxo thoroughly inspects and re-certifies all of its reconditioned biomedical equipment. Batteries are checked or replaced, and each unit is cosmetically repaired. Calibrations are made and all safety inspections will be completed. In most cases our equipment warranties meets or exceeds those of the OEM.