



Edan F6 & F6 Express Fetal Monitor- Fully Refurbished

F6 is a high end hospital grade monitor designed for multi-bed units with central station, yet simple enough for private practice. It offers the most advanced integrated monitoring for the delivery room in big hospital; private obstetrician's office; antepartum clinic. F6 are designed for use in moving situation covers the entire continuum of antepartum and intrapartum applications.

Features:

- Basic parameters: FHR, TOCO, Event Mark, AFM
- Optional Dual FHR monitoring
- Optional Internal parameters: IUP/DECG
- 24 hours waveforms playback
- Optional Built-in wireless module
- 150/152mm wide paper printing
- Support American/International standard
- 1、 2、 3cm/min real-time printing speeds
- Quick printing for stored waveform
- Lithium battery for 4 hours continuous working
- 10.2 - inch high resolution color folding up LCD screen
- Probe rack and wall mounting rack
- Handle for easy carry
- Optional Fetal Stimulators
- Multi-crystal pulse wave Doppler waterproof transducer for FHR detection
- Signals Overlap Verification

Specs:

- Dimensions: 347mm x 330mm x 126mm
- Weight: 6 kg. approx.
- Resolution: 800 x 480
- 14.8 V rechargeable Lithium-ion battery
- Battery Continuous working time: 2-4 hours



Why purchase Auxo Medical Reconditioned Biomedical Equipment?

Auxo thoroughly inspects and re-certifies all of its reconditioned biomedical equipment. Batteries are checked or replaced, and each unit is cosmetically repaired. Calibrations are made and all safety inspections will be completed. In most cases our equipment warranties meets or exceeds those of the OEM.