



Edan F9 & F9 Express Fetal Monitor - Fully Refurbished

F9 is high-end fetal & maternal monitor, it offers the most advanced integrated monitoring for the delivery room in big hospital; private obstetrician's office; antepartum clinic. F9, designed for moving situations, covers the entire continuum of antepartum, intrapartum and postpartum applications.

Features:

- Basic parameters: FHR, TOCO, Event Mark, AFM
- Internal parameters: IUP, DECG (optional)
- 24 hours waveforms playback
- Optional Built-in wireless module
- Software for data transmission to PC
- 150/152mm wide paper
- Support US/International standard paper
- 1,2,3cm/min real-time printing speeds
- Quick printing for stored waveform
- Lithium battery for 4 hours continuous working
- 12.1- inch high-resolution color folding up LCD touch-screen
- Probe rack and wall mounting rack
- Handle for easy carry
- Optional Fetal Stimulator
- 12 - crystal pulse wave Doppler water-proof transducer for FHR detection
- Signal Overlap Verification
- Maternal parameters: MECG, NIBP, MSpO2 are available



Specs:

- Dimensions: 347mm x 330mm x 126mm
- Weight: 6 kg. approx.
- 14.8 V rechargeable Lithium-ion battery
- Continuous Working Time: 4 Hours
- Resolution: 800 (H) x 600 (H)

Why purchase Auxo Medical Reconditioned Biomedical Equipment?

Auxo thoroughly inspects and re-certifies all of its reconditioned biomedical equipment. Batteries are checked or replaced, and each unit is cosmetically repaired. Calibrations are made and all safety inspections will be completed. In most cases our equipment warranties meets or exceeds those of the OEM.