

Edan IM60 - Fully Refurbished

With selectable optional parameters including IBP, cardiac output and capnography, iM60 can adapt to a wide range of sub-acute divisions. The high resolution color TFT-LCD together with the up-to-date customizable shortcut menu design brings easy one-touch access to frequent used functions, allowing more attentions to the patients. Employing low-power consumption design, iM60 is with fanless structure which creates no noise and no dust accumulations.

Being embedded with various communication capabilities including HL7, LAN, built-in Wi-Fi, Nurse Call, and Defibrillator Synchronization, iM60 can easily communicate and co-work with hospital information system and other hospital facilities.

Features:

- 10.4" high resolution TFT-LCD touch screen
- 3/5-lead ECG, NIBP, SpO2, 2-TEMP, PR
- 8.4" color TFT display with maximum 11 waveforms
- Pacemaker detection, electrosurgical interference proof
- Defibrillation protection and defibrillation synchronization
- OxyCRG available to judge the respiration and circulation function for neonates
- Large font display
- Arrhythmia and S-T segment analysis
- Built-in rechargeable Li-ion battery
- Nurse Call function and bi-directional communication with MFM-CMS central station
- HL7 compatible vis XML files
- SD card slot enables memory extension for large data storage
- USB and serial ports enables abundant future upgrades
- VGA output and analog output
- Complete and flexible mounting solution for different scenarios
- Optional Configurations: Thermal Recorder, Touch Screen, Built-in Wi-Fi and CO2 (Sidestream/Mainstream), 2-IBP, C.O.



Why purchase Auxo Medical Reconditioned Biomedical Equipment?

Auxo thoroughly inspects and re-certifies all of its reconditioned biomedical equipment. Batteries are checked or replaced, and each unit is cosmetically repaired. Calibrations are made and all safety inspections will be completed. In most cases our equipment warranties meets or exceeds those of the OEM.