

Edan IM70 - Fully Refurbished


Embedded with the latest improved designs, iM70 is meant to add splendors to the i Series portfolio and answer for the latest clinical requirements across the applications from general cares to anesthesia monitoring, while keeping connected to hospital systems in all divisions. The high resolution color TFT-LCD together with the up-to-date customizable shortcut menu design brings easy one-touch access to frequent used functions, allowing more attentions to the patients. Without making any noise or causing any dust deposition, the no-fan design of iM70 is especially suitable for night time operations and the divisions which require a dust-free environment.

Features:

- 12.1" Color TFT LCD Monitor (Part Number: 03.60.360145)
- 3/5-lead ECG, RESP, NIBP, Edan SpO2, 2-TEMP, PR
- Arrhythmia
- S-T Analysis
- OxyCRG
- Large Font Display
- USB data storage and review
- Bi-directional Communications with Central Monitoring System
- Waveform color adjustable
- Customizable alarm setting

Accessories:

- Skin Temperature Probe
- EDAN SH1 Adult Reusable SpO2 Sensor
- Adult Cuff (25cm-35cm)
- NIBP Tube (3m) with connector
- 3-lead ECG integrative Cable with Leadwires, Snap (AHA) (Defibrillation) (5 lead if indicated)
- Power Cord (USA standard)
- x30 Adult Disposable Adhesive Electrodes (Nonwovens), Snap(30pcs/pack)
- Rechargeable Lithium-Ion Battery (14.8V, 4200mAh)
- User Manual
- Ground Cable



Why purchase Auxo Medical Reconditioned Biomedical Equipment?

Auxo thoroughly inspects and re-certifies all of its reconditioned biomedical equipment. Batteries are checked or replaced, and each unit is cosmetically repaired. Calibrations are made and all safety inspections will be completed. In most cases our equipment warranties meets or exceeds those of the OEM.