



## Edan SE-12 Express ECG - Fully Refurbished

The twelve-channel ECG SE-12 Express is a special system of smart ECG series, which support both resting ECG and Exercise ECG. Alphanumeric keyboard and one touch operation design makes it the perfect choice for you.

### Features:

- 12.1 inch high-resolution foldable color TFT touch-screen
- Specific function keys for Exercise ECG test
- Various treadmills or bicycles and Suntech tango+ exercise blood pressure monitor supported
- Customized exercise protocols
- ST analysis and reports
- Complete digital filters, resisting baseline drift, AC and EMG interference
- Ethernet and RS-232 interface for data transmission to PC
- Built-in rechargeable Lithium battery
- Sleep mode for saving energy and extending LCD life
- Real-time waveform freezing
- 120 seconds ECG waveform reviewable
- Waveforms and auto-diagnosis results preview before printing
- Automatic arrhythmia detection and recording
- 200 ECGs in internal memory
- Automatic measurement and interpretation tested with authoritative CSE database
- Heart rate variability (HRV) analysis
- Alphanumeric keyboard
- Dimensions: 420mm x 330mm x 120mm
- Recorder: Thermal printer
- Record Paper: Folded paper, Rolled paper
- Paper Width: 216mm , 210mm



### Why purchase Auxo Medical Reconditioned Biomedical Equipment?

Auxo thoroughly inspects and re-certifies all of its reconditioned biomedical equipment. Batteries are checked or replaced, and each unit is cosmetically repaired. Calibrations are made and all safety inspections will be completed. In most cases our equipment warranties meets or exceeds those of the OEM.