



Edan SonoTrax Series Fetal Dopplers- Fully Refurbished

SonoTrax Series Ultrasonic Pocket Dopplers provide care givers with the most cost effective solution to meet the routine obstetric and vascular requirements. The SonoTrax Fetal Dopplers offer an easy to use, compact design. The devices also feature a built-in audio recorder, automatic power-off to save power and a large backlight LCD display. Along with these features, the devices have accurate FHR detection with high-fidelity sound effects and are able to support 2, 3, 4, 5 and 8MHz waterproof probes.

Features:

• SonoTrax Lite

- Main Unit (without LCD display)
- 3MHz Waterproof Interchangeable Probe
- Normal Carry Case

• SonoTrax Basic A

- Same as SonoTrax Lite plus:
 - LCD display
 - Blue backlight, battery indication, probe frequency detection
 - 3 Working modes (real-time, average, manual)

• SonoTrax Pro

- Same as SonoTrax Basic A plus:
 - Audio recording and playback

• SonoTrax II

- LCD display
- Blue backlight, battery indication, probe frequency detection
- 3MHz waterproof interchangeable probe
- AC recharger with rechargeable lithium polymer battery (900mAh)
- 3 working modes (real-time, average, manual)
- Normal carry case

• SonoTrax II Pro

- Same as SonoTrax II plus:
 - Audio recording and playback



Why purchase Auxo Medical Reconditioned Biomedical Equipment?

Auxo thoroughly inspects and re-certifies all of its reconditioned biomedical equipment. Batteries are checked or replaced, and each unit is cosmetically repaired. Calibrations are made and all safety inspections will be completed. In most cases our equipment warranties meets or exceeds those of the OEM.