



## GE Aestiva Compact w/ 7100 Ventilator - Fully Refurbished

The GE Aestiva/5 Compact is an advanced anesthesia delivery system with basic ventilation at a moderate price. While maintaining the industry leading excellence of the Aestiva®/5 Anesthesia System - open architecture, integrated breathing circuit, compact design - it features a configurable version of Datex-Ohmeda's newest ventilation platform, the 7100 ventilator. The GE Aestiva 5 offers a unique combination of quality performance and cost-control. The 7100 ventilator provides a wide range of capabilities, including optimal tidal volume compensation and integrated volume and oxygen monitoring.

### Features:

- Smartvent 7100 ventilator
- Patient breathing system with circle module
- Low flow anesthesia
- 3-gas delivery system
- Open, flexible architecture
- Fold up writing shelf
- Lockable drawer
- Light strip
- Integrated Suction
- Integrated O2 flowmeter
- Eight machine hoses/cables integrated
- User configurable drawers/shelving

### Specs:

- Height: 130.3 cm/51.3 in
- Width: 75 cm/29.5 in
- Depth: 83 cm/32.7 in
- Weight: Approximately 106.6 kg/235 lb  
(machine configured with one drawer and no vaporizers)
- 120 x 90 LCD screen



### Why purchase Auxo Medical Reconditioned Biomedical Equipment?

Auxo thoroughly inspects and re-certifies all of its reconditioned biomedical equipment. Batteries are checked or replaced, and each unit is cosmetically repaired. Calibrations are made and all safety inspections will be completed. In most cases our equipment warranties meets or exceeds those of the OEM.