



auxoMEDICAL

GE B40 Patient Monitor - Fully Refurbished

The GE Procure B40 monitor is designed with advanced clinical measurement technologies, including proven GE algorithms, for accurate and reliable monitoring even in challenging conditions. Ideal for monitoring adult, pediatric and neonatal patients, the B40 is modular and configurable to monitor number of possible parameters including non-invasive blood pressure, invasive blood pressure, end-tidal carbon dioxide, heart/pulse rate, respiration rate, ECG, temperature, SpO2, agent identification, as well as advanced gas and respiratory monitoring parameters.

Features:

- Customizable settings and alarms
- Small footprint
- Network connectivity
- 12.1" color display and one-button access to commonly used features.
- Complete portable monitoring solution for adults, pediatric and neonatal patients.
- Graphical 72-hour trend display with numerical data.
- Modular and configurable with advanced parameters.
- Electronic medical record (EMR) system through CARESCAPE server.
- Alarm system with preset limits for critical alerts.
- Easy to use menus to save on time and training.
- Multiple mounting options
- Rugged design stands up to heavy use in busy care areas



Why purchase Auxo Medical Reconditioned Biomedical Equipment?

Auxo thoroughly inspects and re-certifies all of its reconditioned biomedical equipment. Batteries are checked or replaced, and each unit is cosmetically repaired. Calibrations are made and all safety inspections will be completed. In most cases our equipment warranties meets or exceeds those of the OEM.