



## Hausted 493 - Fully Refurbished

The Hausted Fluoro-Track Fluoroscopy Capable Stretcher is ideal for fluoroscopy, endoscopy, and other procedures. The Fluoro-Track has a C-Arm clearance of 24" and has a patient weight capacity of 500 lbs. The 493 features Sliding Top that travels a total of 12" for positioning during procedures. Radiolucent Patient Surface. It travels easily on 8" Casters with Four Brake/Steer Levers located at each corner for easy implementation and steers from all four corners

### Features:

- Stainless Steel 6-Leg, 3-Position Retractable Rails with Built-In Push-Steer Handles and Top Caps that Eliminate
- Patient Transfer Gap
- Dual Side Foot Pedal Hydraulic Controls for Height Adjustment 23° - 36°, Trendelenburgs: 0° - 15°
- Pneumatic Assist Backrest Adjusts 0° - 80°
- Brake and Steering Activation from All 4 Corners
- 8" Tente Casters
- Powder Coated Steel Patient Surface
- Seamless Thermoformed Base Cover with Storage and Built-In Oxygen Tank Holder
- Full Perimeter Non-Marring Bumper System
- Full-Length Radiolucent Top Meets BRH Requirements
- Translating Top with 18" Movement 12" Towards Foot and 6" Towards Head for Patient Positioning during
- Fluoroscopic Procedures
- Patient Surface: Nominal 24" x 76"
- Weight Capacity: 500 lbs.



### Why an Auxo Medical Refurbished Table or Stretcher?

Auxo tables and stretchers are refurbished to the some of the highest standards in the industry.

We completely disassemble each table and stretcher, and completely rework or replace all worn parts and return each table and stretcher to OEM standards. Each table and stretcher is reassembled, tested and is warranted against defects.