



Hill-Rom Centra Series 850 & 852- Fully Refurbished

Engineered to meet the highest safety standard for electric medical equipment, the bed combines aesthetically pleasing designs with innovative features to create an ideal bed for any long-term care environment. The Hill-Rom® Centra® Series 850 & 852 hospital beds are perfect for hospitals & surgery centers. These beds offers a variety of features that are perfect for your patient care needs.

Features:

- Three motors
- High-low travel time: 21 seconds
- Sleep surface: 36" x 82"
- Bed frame clearance with 5" casters: 7.5"
- Low position: 17"
- High position: 30"
- Central pedal for braking and steering
- Night light
- CPR Releas
- Retractable Frame
- Parallelogram Side rails (optional)
- Controls on the foot end and on the side rails
- Electrically adjust height, knee & back sections
- Foot section can be raised manually
- Side rail controls can be locked out by care giver
- Separate head and foot end side rails
- Trendelenburg position to facilitate ease of care and add comfort (the entire bed is placed 12 degree angle)



Why an Auxo Medical Refurbished Table or Stretcher?

Auxo tables and stretchers are refurbished to the some of the highest standards in the industry.

We completely disassemble each table and stretcher, and completely rework or replace all worn parts and return each table and stretcher to OEM standards. Each table and stretcher is reassembled, tested and is warranted against defects.