



Lifepak 9 With Pacing - Fully Refurbished

The Lifepak 9 from Physio-Control, a microprocessor-based unit, is available in models: 9, 9P (pacing), and 9A. These units have a built-in recorder and sealed lead acid battery back-up. Lifepak 9s meet current ACLS standards and allows for 17 energy selections, HR alarms, 5-lead monitoring, recorder and CODE SUMMARY critical event record. The Lifepak 9 is defibrillator/monitor, cardioscope and non-invasive pacemaker. It provides 17 energy selections, heart rate alarms, three-lead ECG monitoring, recorder and Code Summary critical event recorder.

Features:

- Internal battery supply
- Available with optional pediatric paddles, internal paddles, pacing cable and pads, and ECG patient cable and leads.
- Isolated ECG via Quick-Look defibrillator paddles, Fast-Patch disposable defibrillation/ECG electrodes or three-lead patient cable.
- Operates on either AC Power or Internal Battery power
- Sealed lead acid battery back-up
- Built-in recorder
- Pacing
- Meet current ACLS standards and has annotating summary
- Allows for 17 energy selections, HR alarms, 5-lead monitoring, recorder and CODE Summary critical event record
- Provides:
 - Defibrillation
 - Synchronous Cardioversion
 - Cathode Ray Tube monitor screen
 - A strip chart recorder
 - 3-lead electrocardiogram monitoring
 - Heart rate display and alarms
 - Multiple energy selection
- Dimensions: 35.2 cm(H) x 29.7 cm(W) x 31 cm(D), 13.9" x 11.7" x 12.2"
- Weight: 12.7 kg (28 lbs)



Why purchase Auxo Medical Reconditioned Biomedical Equipment?

Auxo thoroughly inspects and re-certifies all of its reconditioned biomedical equipment. Batteries are checked or replaced, and each unit is cosmetically repaired. Calibrations are made and all safety inspections will be completed. In most cases our equipment warranties meets or exceeds those of the OEM.