



# auxoMEDICAL

## Medivator SSD-102LT- Fully Refurbished

Automated endoscope reprocessor features one reprocessing basin with user-programmable reprocessing cycles, and an enhanced sensor package for continuous monitoring of all cycle parameters. The SSD-102LT includes automated endoscope leak testing, a user-selectable wash cycle for endoscope cleaning, automatic alcohol purge and air drying of endoscope channels, and an optional data management system for electronic cycle tracking and record keeping

### Features:

- Fast cycle time for quick endoscope turnaround
- Cycle parameters continuously monitored
- Leak test ensures endoscope protection throughout cycle
- Automated alcohol flush and air dry
- Audible and visual alarms for reprocessing cycle deficiency
- Assures patient and operator safety
- Cycle customization meets specific reprocessing needs
- Bacteria retentive three-stage water filtration system
- Compatible with all current models of endoscopes, including those with elevator wire channels
- Printed documentation of all reprocessing cycle data
- Self-disinfects internal fluid lines during every cycle
- For use with most FDA-cleared reusable high-level disinfectants
- Computerized endoscope reprocessing records, vapor management system and neutralization system optional

### Specs:

- Dimensions: 46" H x 22" W x 21" D
- Height with lid open: 65"(165 cm)
- Weight: 220 lbs (100 kg)
- Reservoir Capacities: Disinfectant: 4 gallons (15.1 liters) Alcohol and detergent: 27oz each (800 ml)
- Disinfectant Temperature Range: Room temp to 125°F (52°C)



### Why purchase Auxo Medical Reconditioned Biomedical Equipment?

Auxo thoroughly inspects and re-certifies all of its reconditioned biomedical equipment. Batteries are checked or replaced, and each unit is cosmetically repaired. Calibrations are made and all safety inspections will be completed. In most cases our equipment warranties meets or exceeds those of the OEM.