



Midmark 404 Exam Table- Fully Refurbished

The Midmark® 404 Exam Table is an efficient and comfortable exam table designed with the patient and physician in mind. This table has an easily adjustable top for positions from horizontal to 82 degrees, a leg extensions with concealed pan, a removable upholstered top for easy cleaning, lots of storage and infinitely adjustable stirrups. This table is sold as refurbished in excellent condition and is available in a variety of colors.

Features:

- Height in table position: 31.75"
- Table length: 54" extends to 72"
- Table Width: 27"
- Foot extension length: 18"
- Foot extension width: 20"
- Stirrups
- Pneumatic Back
- Storage Drawers
- Contoured foam construction with antimicrobial vinyl cover
- Drawers (3 Side/3 Front)
- Pull out step on Front
- Max Patient Weight 300lbs
- Back Section Adjustment: 0° to 80°
- Pelvic Lift Adjustment: Seat Section Adjustment to +5°
- Stirrup Extension (Max): 14in
- Stirrup Lateral Positions:
- (Max): 20.5in
- Footrest Extension Adjustment: to 7°
- Dimensions: 27in W x 55in L x
- Weight: 290 lbs
- Paper Rolls Holder: 21.0" x 3.5" rolls (53.3 cm x 8.9 cm).
- Treatment Pan (Polystyrene): 9.5" x 9.5" x 3.0" (24.1 cm x 24.1 cm x 7.6 cm).
- Electrical Rating: 115 VAC, 60 Hz, 5 Amps.



Why purchase Auxo Medical Reconditioned Biomedical Equipment?

Auxo thoroughly inspects and re-certifies all of its reconditioned biomedical equipment. Batteries are checked or replaced, and each unit is cosmetically repaired. Calibrations are made and all safety inspections will be completed. In most cases our equipment warranties meets or exceeds those of the OEM.