



## Midmark 417 Podiatry Chair - Fully Refurbished

The 417 has everything you'd want in a chair custom-built for your practice: contemporary styling with clean, handsome lines – and feature-packed convenience for podiatrist, patient and staff! Its foot control couldn't be easier to program, or more responsive. A touch is all it takes to position each patient's feet precisely where you want them. Or re-position them during procedures without compromising your sterile field. And it won't slip, slide – or scuff the floor. The optional Hand-held Commander mounts on either side of the foot section to provide 7 additional choices.

### Features:

- Adjustable height, tilt, – work standing or sitting
- Fastest power podiatric chair – 1 button control
- Foot control (standard) – 3 power positioning, 2 prog choices
- One-Touch automatic return – saves seconds on each procedure
- Deep, sculptured cushions – built-in pillow
- Duplex hospital-grade electrical receptacle – left and right
- Vinyl upholstery – stain-resistant; Velcro-fastened
- ISO 9001 certified
- Powder epoxy-baked paint finish; Black base
- Painted steel surface/molded housings: Pearl Grey
- Height, tilt and back move simultaneously
- Arm rests rotate 90¼ for easy chair entry/exit



### Specs:

- Convenient low height (22.5")
- High range of 41" (the foot section can be as high as 58" when the table is at maximum height and in tilt)
- Width (Back Section) - 22.75" (57.8cm)
- Width (Seat Section) - 22.75" (57.8cm)
- Width (Foot Section) - 21.25" (53.9cm) tapering to 17.75" (44.4cm)
- Electrical - Domestic: 115 VAC, 60Hz, 12AMP

### Why purchase Auxo Medical Reconditioned Biomedical Equipment?

Auxo thoroughly inspects and re-certifies all of its reconditioned biomedical equipment. Batteries are checked or replaced, and each unit is cosmetically repaired. Calibrations are made and all safety inspections will be completed. In most cases our equipment warranties meets or exceeds those of the OEM.