



Midmark Ritter M11 - Fully Refurbished

The chamber size gives the Ritter M11 the distinction of having the largest sterilizing capacity of any standard countertop sterilizer on the market, making it perfect for larger packs and instruments. Ritter sterilizers are designed to be easier and are time tested for reliability. With any Ritter sterilizer, simply load it, set it and forget it. Select the cycle and the time is set for you. Press "Start" and the sterilizer will fill with the correct amount of water. Once the sterilization process is complete, the door opens automatically and quietly to dissipate the steam, giving you a fast drying time for your instruments.

Features:

- Exclusive automatic-opening door for ease of use and a quick drying time
- Large 11" by 18" chamber for increased capacity
- Steam-Flush Pressure Pulse Air Removal System ensures proper steam penetration
- Self Program Controls enable you to create different cycle parameters to meet specific instrument processing needs
- Pre-Programmed Controls for convenience
- Green LCD Display for easy view of indicators
- Configured for plug-in or hard-wired connection
- All sterilization parameters are digitally displayed



Specs:

- Overall length w/plug: 23.8" (60.5 cm)
- Overall width: 17.8" (45.2 cm)
- Overall height w/printer: 17.8" (45.2 cm)
- Minimum Countertop area: 17.8" x 21.0" (45.2 x 53.3 cm)
- (2) Large Tray Size - 9" W (22.9 cm) x 15" L (38 cm) x 1.1" D (2.9 cm)
- (2) Small Tray Size - 6.6" W (16.8 cm) x 15" L (38 cm) x 1.1" D (2.9 cm)
- Unit Weight: 99 lbs. (44.9 kg)
- 11" x 18" chamber

Why purchase Auxo Medical Reconditioned Biomedical Equipment?

Auxo thoroughly inspects and re-certifies all of its reconditioned biomedical equipment. Batteries are checked or replaced, and each unit is cosmetically repaired. Calibrations are made and all safety inspections will be completed. In most cases our equipment warranties meets or exceeds those of the OEM.