



## Midmark Ritter M9 - Fully Refurbished

The Ritter M9 is an autoclave that couldn't be easier and is time tested for reliability. You load it set it and forget it. Select your cycle and the time is set for you. Press Start and the Midmark M9 fills with the correct amount of water. Once the sterilization process is complete, the door opens automatically and quietly to dissipate the steam to give you fast drying time for your instruments. Its real convenience!

### Features:

- Fills with correct amount of water when you hit the start button.
- Door opens automatically to dissipate steam and give you the fastest drying time.
- Built-in safety features such as overheat protection, a pressure relief valve and flashing "Error" light.
- Self Program Controls enable you to create different cycle parameters to meet specific instrument processing needs
- Green LCD Display for easy view of indicators
- Overall Dimensions: 18.81"(L) x 15"(W) x 16.44"(H)
- Countertop area: Minimum depth 18" (45.7 cm)
- Chamber: 9" diameter (22.9 cm) 15" depth (38.1 cm)  
3.5 gal. usable volume (13.4 liters)
- Trays: Two large - 12" L x 7.31" W x .88" D (30.5 x 18.6 x 2.2 cm) Two small - 12" L x 5.62" W x .88" D (30.5 x 14.3 x 2.2 cm)
- Unit Weight: 70 lbs. (31.8 kg)
- Water reservoir capacity: .88 gal. (3.3 liters)
- Electrical 120 VAC, 60 Hz, 12 AMPS 120 VAC has 1425 watt tubular immersion heater
- Unwrapped - 14 min. Hot, 19 min. Cold, 3 min. Sterilization Time
- Pouches - 17 min. Hot, 22 min. Cold, 5 min. Sterilization Time
- Packs - 40 min. Hot, 46 min. Cold, 30 min. Sterilization Time
- Open Door Dry Cycle



### Why purchase Auxo Medical Reconditioned Biomedical Equipment?

Auxo thoroughly inspects and re-certifies all of its reconditioned biomedical equipment. Batteries are checked or replaced, and each unit is cosmetically repaired. Calibrations are made and all safety inspections will be completed. In most cases our equipment warranties meets or exceeds those of the OEM.