



Mindray A5 Anesthesia Machine - Fully Refurbished

The Mindray A5 Anesthesia Machine provides today's clinicians with a variety of advanced ventilation modes, an easy-to-use interface and an ergonomic design. The A5 features a 15" touch screen that allows operators to navigate easily between ventilation modes, allowing caregivers to spend more time with patients. The A5 is capable of outputting data in HL7 protocol, the industry-recognized format for EMR systems.

Features:

- Volume Control Ventilation (VCV)
- Pressure Control Ventilation (PCV) with volume guarantee
- Pressure Support Ventilation (PS)
- Synchronized Intermittent Mandatory Ventilation in volume and pressure modes (SIMV-VC and SIMV-PC)
- Automatic compliance compensation and fresh gas compensation maintain accurate tidal volumes
- The warmed Breathing System virtually eliminates internal condensation
- The built in caster guard design clears away cables and hoses
- The economical single container absorber reduces compressible dead space volume and accepts non-proprietary prepackaged and loose fill absorbent
- Integrated backup control touchpad
- Deck lighting with adjustable brightness control
- Two hour battery backup
- Central brake controls all 4 casters

Specs:

- Overall Dimensions: 55.1"(H) x 41.3" (W) x 31.7"(D)
- Weight (without vaporizers and gas cylinders) 353 lbs (160 kg)
- Top Shelf Weight Limit: 88 lbs (40 kg)
- Top Shelf Dimensions: 24.2"(W) x 14.2"(D)
- Work Surface Dimensions: 24.2"(W) x 15"(D) x 33.4"(H)
- Drawer Dimensions: 5.3"(H) x 17.3"(W) x 15.1"(D)



Why purchase Auxo Medical Reconditioned Biomedical Equipment?

Auxo thoroughly inspects and re-certifies all of its reconditioned biomedical equipment. Batteries are checked or replaced, and each unit is cosmetically repaired. Calibrations are made and all safety inspections will be completed. In most cases our equipment warranties meets or exceeds those of the OEM.