



NUVO Surgical V800 - Fully Refurbished

Designed to accommodate surgical procedures in nearly all fields, including urology, neurology, orthopedics, and cardiology, the Nuvo V800 series surgical tables set new standards for versatility, innovative features and solid dependability. With design and engineering thought through to a remarkable degree, they provide such benefits as unprecedented C-arm access, self-diagnosis with Instant Error Code display, and a control system that automatically prevents of sectional conflicts.

Features:

- 800 lb. Weight Capacity (V800)
- Removable, Split Leg Section (V800)
- Reverse Orientation
- Slide and Hybrid Options
- Override Controls
- Perineal Cutout
- Table Leveling
- Height: 26.8 - 44 in (68 - 112 cm)
- Lateral Tilt: 20°
- Trendelenburg and Reverse Trendelenburg: 30°
- Back Section: 80° up or 50° down
- Leg Section: 30° up or 90° down
- 180° open to the left and right.
- Head Section: 90° up or 45° down
- Kidney Elevator Elevation 0 ~ 4.7 in (12cm)
- Flex: Maximum of 25° down
- Reflex: Maximum of 30° up

Also Available with the following features:

- Kidney Bridge
- Floor Locks and Kidney Bridge
- Floor Locks, Sliding Top, and Kidney Bridge
- Kidney Bridge and Sliding Top



Why an Auxo Medical Refurbished Table or Stretcher?

Auxo tables and stretchers are refurbished to the some of the highest standards in the industry.

We completely disassemble each table and stretcher, and completely rework or replace all worn parts and return each table and stretcher to OEM standards. Each table and stretcher is reassembled, tested and is warranted against defects.