



Narkomed Mobile Anesthesia Machine- Fully Refurbished

The Drager Narkomed M is a state-of-the-art anesthesia machine that is ideal for efficient transport and quick assembly, making it ideal for military use and other mobile environments that are tight on space. The anesthesia machine's design is simple with fewer components than industry competitors, allowing medical personnel to assemble and disassemble the Narkomed M in labor and time saving manner. The modular architecture is durable and versatile.

Features:

- System design minimizes the number of required components thereby minimizing assembly/disassembly time and reducing loss of components.
- Modular architecture is extremely robust and flexible enough for use in a wide variety of hospital settings.
- Stored and transported in its own custom air tight containers.
- Simple, modular design is easy to assemble
- Small footprint ideal in tight spaces
- Quick disassembly for easy transport
- Closed (Low Flow) anesthesia techniques
- Semi-closed anesthesia techniques
- Opened anesthesia techniques
- Break-Resistant Flow Meters
- Integrated Ventilator: flow rate of 0 - 5 LPM
- Adaptable to use in both active and passive gas scavenging systems



Specs:

- Dimensions: 23.6" (59.9 cm) W x 24.4" (62.9 cm) L x 53" (134.6 cm) H
- Weight: 75 lbs (34 kg)
- Power: 110 / 220V, 50/60Hz, 250 watts, 2.5 amps max
- Battery back-up to fully power all operations for a minimum of one hour (brand new battery installed)

Why purchase Auxo Medical Reconditioned Biomedical Equipment?

Auxo thoroughly inspects and re-certifies all of its reconditioned biomedical equipment. Batteries are checked or replaced, and each unit is cosmetically repaired. Calibrations are made and all safety inspections will be completed. In most cases our equipment warranties meets or exceeds those of the OEM.