



auxoMEDICAL

OEC 6600 Mini C-Arm - Fully Refurbished

The OEC 6600 Mini C-Arm brings the quality and reputation you have come to expect from the world leader in C-Arm design. The 6600 is a portable C-Arm system that can be used for the imaging of extremities. It comes with the option of a 4" or 6" image intensifier and two 16" monitors. It is easy to use and small enough to be operated by a staff person or surgical technician. It does not need to be operated by an X-ray technologist. Depending on its configuration, the OEC 6600 Mini C-Arm can store a minimum of 64 to a max of 200 images.

Features:

- 4/6" Tri Mode Image Intensifier
- Rotating Anode X-Ray Tube
- Digital Camera
- GSP – General Surgical Package Software
- LIH – Last Image Hold
- Boost Function
- Sharpen Feature
- High Frequency Generator
- Set To Factory OEC Specs
- Sony Paper Printer
- Patient Annotation Keyboard
- Dual 16 Inch Hi-Resolution Monitors 6 Inch Image Intensifier
- Surgery Platform
- Image Storage
- Motion Artifact Reduction System Gamma Correction
- One Shot Frame Integration
- 200 image storage Frame averaging Patient annotation Auto Fluoro Mode Keyboard
- Reverse camera keys Footswitch
- 27" x 27" footprint



Why purchase Auxo Medical Reconditioned Biomedical Equipment?

Auxo thoroughly inspects and re-certifies all of its reconditioned biomedical equipment. Batteries are checked or replaced, and each unit is cosmetically repaired. Calibrations are made and all safety inspections will be completed. In most cases our equipment warranties meets or exceeds those of the OEM.