



auxoMEDICAL

OEC 9600 C-Arm- Fully Refurbished

A platform from which you can expand and add capabilities as desired, the GE OEC 9600 mobile c-arm is an imaging tool that is capable of performing different roles in multiple environments. Light and compact for easy transportation, the OEC 9600 nevertheless has a large depth of 26 inches and a more-than-usual amount of space between the x-ray tube and image intensifier. The extra space allows for more room to perform procedures and to accommodate protruding instruments. It can then be used for laparoscopic and endoscopic procedures. The OEC 9600 Digital Mobile C-Arm provides great images with low dose exposure. The 9600 is maneuverable, is simple to operate, & provides large images

Features:

- Rotating X-Ray Tube
- Tri-Mode 9/6/4.5 inch Image Intensifier
- High Resolution CCD
- Fluoro Boost
- Pulsed Fluoro
- Image Storage
- Dual Hi-Resolution Monitors
- Last Image Hold
- Motion Artifact Reduction
- Edge Enhancement
- Sharpen Function
- Gamma Correction
- Floppy Drive
- Dual 16-inch high resolution square monitors
- 4.0 kW generator
- Workstation Dimensions: 27.25"(W) x 68.25"(H) x 27.25(D)
- C-arm Dimensions: 78.5"(L) x 68.25"(H) x 33"(W)



Why purchase Auxo Medical Reconditioned Biomedical Equipment?

Auxo thoroughly inspects and re-certifies all of its reconditioned biomedical equipment. Batteries are checked or replaced, and each unit is cosmetically repaired. Calibrations are made and all safety inspections will be completed. In most cases our equipment warranties meets or exceeds those of the OEM.