



Pelton & Crane Magnaclave - Fully Refurbished

The Pelton & Crane Magnaclave is a fully self-contained steam sterilizer with a large chamber, making it ideal for busy group practice, emergency clinics, surgery centers and small hospitals. The chamber capacity allows it to hold up to 14 Hu-Friedy cassettes or equivalent. The unit is sold with on basket tray and one medium shallow tray. It is equipped with a safety valve which limits the maximum chamber pressure to 38 PSI. This unit is sold refurbished and includes a stand. Electrical Requirements: Dedicated 30 amp circuit/220-240 VAC or 208 VAC.

Features:

- Safety lock system protects personnel from accidentally opening the door under pressure.
- Completes a flash cycle in 18 minutes from a warm start, 35 minutes from a cold start. Drying phase not included
- Reservoir drain is conveniently located on the front for easy cleaning.
- Simple and easy door locking mechanism.
- Water reservoir has large capacity to adequately cool condensed steam, reducing the need for refilling.
- Water level indicator
- Preheat feature allows faster sterilization time.
- Optional Stand available
- Timer activates automatically as the desired pressure for sterilization is reached

Specs:

- Large chamber capacity (16" dia x 30" deep)
- Total Dimensions: 28"W, 34.5"D, 23" Tall (53" Tall w/ Stand)
- Weight: 355 lbs (the actual sterilizer), 89lbs (stand only)



Why purchase Auxo Medical Reconditioned Biomedical Equipment?

Auxo thoroughly inspects and re-certifies all of its reconditioned biomedical equipment. Batteries are checked or replaced, and each unit is cosmetically repaired. Calibrations are made and all safety inspections will be completed. In most cases our equipment warranties meets or exceeds those of the OEM.