



Ritter M11 Autoclave, Old Style - Fully Refurbished

Refurbished Ritter M11 with 11" x 18" chamber. It's chamber size gives it the distinction of having the largest sterilizing capacity of any standard countertop sterilizer on the market, making it perfect for larger packs and instruments. This product is primarily used by veterinarians or surgeons, but is also used in a range of other settings, including active surgery center, medical practice, dental office, schools, industrial lab, doctor's clinic, or tattoo establishment.

Features:

- Exclusive automatic-opening door for ease of use and a quick drying time
- Large 11" by 18" chamber for increased capacity
- Steam-Flush Pressure Pulse Air Removal System ensures proper steam penetration
- Pre-Programmed Controls for convenience
- Red LCD Display for easy view of indicators
- Overall Length with Plug: 23.8" (60.5 cm).
- Overall Width: 17.8" (45.2 cm).
- Overall Height with Printer: 17.8" (45.2 cm).
- Minimum Countertop Area: 17.8" x 21.0" (45.2 cm x 53.3 cm).
- Chamber:
 - 11" Diameter (28 cm).
 - 18" Depth (45.7 cm).
 - 6.5 gallons usable volume (24.6 liters).
- Trays:
 - 2 Large: 9.0" (W) x 15.0" (L) x 1.1" (D) (22.9 cm x 38.0 cm x 2.9 cm).
 - 2 Small: 6.6" (W) x 15.0" (L) x 1.1" (D) (16.8 cm x 38.0 cm x 2.9 cm).
- Unit Weight: 99 lbs (44.9 kg).
- Shipping Weight: 131 lbs (59.4 kg).
- Water Reservoir Capacity: 1.4 gallons (5.3 liters).
- Electrical Requirements:
 - 115 VAC, 50/60 Hz, 15 amps, single phase.
 - 115 VAC has 1425 watt tubular immersion heater.



Why purchase Auxo Medical Reconditioned Biomedical Equipment?

Auxo thoroughly inspects and re-certifies all of its reconditioned biomedical equipment. Batteries are checked or replaced, and each unit is cosmetically repaired. Calibrations are made and all safety inspections will be completed. In most cases our equipment warranties meets or exceeds those of the OEM.