



Scican Statim 5000 G4 Cassette - Fully Refurbished

The fully automatic, Statim 5000 G4 utilizes four, single touch sterilization programs to optimize instrument cleaning, from start to finish. The Statim 5000 G4 has an additional sterilization cycle and a roomier cassette in order to process larger instruments. It's reliable and compact enough to fit in patient treatment areas and fast enough to sterilize expensive Lasik, PRK, endoscopic and cataract instruments, including phaco handpieces, between procedures. In just 17.5 minutes the Statim 5000 G4 can complete a wrapped cycle and only requires 9 minutes to complete an unwrapped sterilization cycle. The Statim 5000 G4 series has been modernized with a new color touch screen and intuitive controls.

Features:

- A large 3.5" high-resolution touchscreen offers a vivid display of messages and current cycle information all with extraordinary clarity.
- Ultra-fast 6 minute unwrapped cycle
- 14 minute wrapped cycle
- Exceptional instrument turnover
- 10 minute Dri-Tec Drying Technology
- Aseptic transfer of instrument directly to the point-of-use
- Reduced instrument investment
- Prolongs instrument life
- Machine Size: 21.5"x16.25"x7.5"
- Cassette Internal Dimensions: 15"x7"x3"
- Reservoir capacity 4 liters
- Weight: 33 kg / 72 lbs
- Power: 11 - 120v 50/60 Hz 1300 W
- Built-in data logger
- Can copy data to the built in USB, view it on the screen or store in on iCloud
- Removable and Portable Cassette



Why purchase Auxo Medical Reconditioned Biomedical Equipment?

Auxo thoroughly inspects and re-certifies all of its reconditioned biomedical equipment. Batteries are checked or replaced, and each unit is cosmetically repaired. Calibrations are made and all safety inspections will be completed. In most cases our equipment warranties meets or exceeds those of the OEM.