



## Single Bay Surgical Sink W/ Knee Operation - Fully Refurbished

This high quality medical scrub sink is great choice for value and price. Fitted with a knee activated control, the single bay scrub sink is a cost effective way to upgrade or outfit your surgical facility. They have a deep basin and a high mount swivel gooseneck spout with adjustable swivel spray aerator for easy and convenient use. The sink is crafted from 304 stainless steel and polished to a #4 finish for a sleek look and easy clean up.

### Features:

- Deep sloping basin
- Adjustable thermostatic controlled pressure regulating mixing valve, one per station
- Check stop strainers for ease of maintenance
- High mount swivel gooseneck spout with adjustable swivel spray aerator
- Flow control valve
- Backflow preventers
- Infrared sinks are equipped with a 24 volt plug in transformer
- Splash screen dividers with double and triple station sinks
- Main support and Z-brackets (in-wall support optional)
- 304 stainless steel polished and ground to a #4 finish



### Why purchase Auxo Medical Reconditioned Biomedical Equipment?

Auxo thoroughly inspects and re-certifies all of its reconditioned biomedical equipment. Batteries are checked or replaced, and each unit is cosmetically repaired. Calibrations are made and all safety inspections will be completed. In most cases our equipment warranties meets or exceeds those of the OEM.