



Skytron 7000 Micro Procedure - Fully Refurbished

Every aspect of the Skytron 7000 Series has been optimized for Microsurgery. The table combines all the traditional degrees of articulation for general surgery with a totally unique 4-way lateral/ longitudinal slide to facilitate microscope procedures. Not only does the table offer 4-way slide, it also features a return-to-center feature that is operated via the custom designed foot control. The Skytron 7000 Micro Surgical Table is designed with a lateral and longitudinal top sliding capability. 2 in any direction eliminates moving, adjusting and refocusing the microscope during micro-surgical procedures.

Features:

- Lateral and longitudinal top slide (2 any direction)
- Two-piece back section with removable upper piece
- Radiotranslucent top throughout
- Removable split leg sections swing out or drop down as needed
- Foot control permits the surgeon to operate slide and return functions for interface with microscope.
- Pendant control of all functions
- Standardized slide rails
- Low height (23") height adjustment
- X-Ray cassette tunnel
- X-Y movement
- Speed adjustment and return



Specs:

Height: 23"-44"

Trendelenberg: 45°

Reverse Trendelenberg: 20°

Weight Limit: 400 Lbs

Dimensions: 76"(L) x 20"(W)

Why an Auxo Medical Refurbished Table or Stretcher?

Auxo tables and stretchers are refurbished to the some of the highest standards in the industry.

We completely disassemble each table and stretcher, and completely rework or replace all worn parts and return each table and stretcher to OEM standards. Each table and stretcher is reassembled, tested and is warranted against defects.