



Skytron Stellar Series Lights - Fully Refurbished

Stellar's high performance and extremely durable polymer Surgical Lightheads produce high intensity, soft white surgical lighting when and where it is needed most, while providing effortless movement. Sealed lightheads and smooth lines protect against dust accumulation, providing quick and easy cleanings between surgical cases.

Features:

- Advanced Polymer Construction
- Bulb life is generally 5X competitive systems, protected by Skytron's Soft-Start System
- CRI of 93 provides tissue differentiation
- Color temperature: 4,000 K
- Featherlight maneuverability
- 360° radial arm and surgical lighthead rotation
- 360° pitch, roll and rotational axis
- 20° upward travel for taller surgeons
- 90° downward travel for low angle lighting
- Multi-light source design protects against surgical interruption during bulb failure
- 23" or 29" Lightheads depending on model
- Camera systems are available with all lighting configurations



23" Models:

- ST23 — Single Radial Arm (130,000 total lux)
- ST2323 — Dual Radial Arms (260,000 total lux)
- ST232323 — Triple Radial Arms (390,000 total lux)

29" Models:

- ST29 — Single Radial Arm (160,000 total lux)
- ST2923 — Dual Radial Arms (290,000 total lux)
- ST292323 — Triple Radial Arms (420,000 total lux)

Why purchase Auxo Medical Refurbished Lights?

Auxo Lights, halogen and LED, are refurbished to the highest standards in the industry. We completely disassemble each light and arm assembly, and completely rebuild or replace all worn parts and return each assembly to OEM standards. Each Light is reassembled, tested and is warranted against defects.