



auxoMEDICAL

Steelco DS 610 - Fully Refurbished

Steelco DS 610 is a high capacity model, frontal loading or pass-through washer disinfector, with HEPA-filtered forced air drying distribution on double main for a perfect distribution on all chamber parts and basket levels. The ergonomic design of the door level height allows the user a convenient loading/unloading job with the support of a manual loading/unloading trolley.

Features:

- High safety and easy to use soft touch control system on full glass panel with LCD color display to visualize cycle status.
- 20 standard pre-programmed cycles available and 20 additional cycles available for customization.
- RS 232 port for printer connection to monitor and validate washing phases..
- Washing and DI disinfection temperature are fully adjustable up to 93°C.
- Heavy washing pump ensures high water flow rate combined with effective spray pressure.
- Washing/drying injection system on the same circuit.
- Inner cabinet, washing arms and tank filters entirely made of high quality AISI 316L stainless steel
- Triple water filtering system
- Two standard automatic liquid dispensers provided with level control.
- Large storage on the basement to allocate up to three 5lt (1.32 Gal US) chemical containers.

Specs:

Dimensions WxDxH: 25.59" x 27.04" x 72.44" (standard configuration single door model)

Chamber volume ~250 lt / 66.04 Gal US

Basket volume ~220 lt / 58.12 Gal US



Why purchase Auxo Medical Reconditioned Biomedical Equipment?

Auxo thoroughly inspects and re-certifies all of its reconditioned biomedical equipment. Batteries are checked or replaced, and each unit is cosmetically repaired. Calibrations are made and all safety inspections will be completed. In most cases our equipment warranties meets or exceeds those of the OEM.