



Steris Amsco M70 Series - Fully Refurbished

The Steris Amsco M70 Series line of blanket warmers was one of the most popular lines ever manufactured. This series features a tabletop version with a rotary knob that controls temperature. The Amsco® M70® also features over-temperature alarms (visual and audible) and a heating indicator (visible). On full size models, the upper compartment of the Amsco® M70® is the same as the tabletop version, but with a separate lower cabinet that is larger in size than the upper.

Features:

- Temperature sets to 120°C
- Door with rubber gasket, /li>
- Two adjustable stainless steel shelves in the bottom compartment
- Easy grip handle mounted at bottom of door
- Heat and overtemp lights
- Five foot the power cord with hospital grade plug
- LED readout for instant temperature confirmation
- Superior durable materials to withstand years of use
- Reliable heating with uniform and efficient heat transfer
- A design to meet AORN Recommended Practices for a Safe Environment of Care, including technology to warm IV fluids (RP VIII.b.3)
- Temperature lock-out for consistency and control
- Celsius and Fahrenheit selection for convenience and accuracy
- Overall Dimensions: 30" Wide x 20.5" Deep x 75.75" High
- Top Chamber Dimensions: 27"x 17" x 14"
- Compartments: 1 o 2
- Installation: Free-standing or recessed
- Visual access to inventory with dual-pane tempered glass or traditional solid steel doors



Why purchase Auxo Medical Reconditioned Biomedical Equipment?

Auxo thoroughly inspects and re-certifies all of its reconditioned biomedical equipment. Batteries are checked or replaced, and each unit is cosmetically repaired. Calibrations are made and all safety inspections will be completed. In most cases our equipment warranties meets or exceeds those of the OEM.