



auxoMEDICAL

Steris Reliance Synergy - Fully Refurbished

The Reliance Synergy Washer/Disinfector is a mechanical washer equipped with an Eagle 3000 Stage 3 control system. The washer/disinfector is designed with nine adjustable cycles. One additional cycle is available for customized programming to meet specific operating requirements. The Steris Reliance® Synergy® Washer / Disinfector is designed to meet increasing reprocessing needs of hospital CSSD's and Ambulatory Surgery centers through improved efficiency and reduced utility consumption.

Features:

- Faster, optimized cycles
- Touch screen, PC controls can interface with instrument tracking and OR scheduling
- Longer spray arms and outward spraying nozzles
- Large capacity
- Improved drying system
- Integrated flow meters provide verification of chemistry injection rates, even down to the Prolystica® Ultra Concentrate levels
- Easy-to-clean spray arms; simple assembly
- Self-flushing centrifugal filter system prevents spray arms from clogging and helps minimize maintenance
- Racks are capable of lumen flushing
- Cool-touch handles on racks for safe handling of racks once they exit the washer
- Environmentally friendly system uses less water, detergent and utilities
- The system comes with one enzyme and one detergent injection pump
- Double-walled, insulated construction
- Pulsed enzyme treatment
- Airtight doors and HEPA filter



Specs:

- Exterior dimensions: 42x80x36-1/2" (1067x2032x927mm)
- Interior chamber dimensions: 26-1/2 x 26-1/4 x 26" (673 x 667 x 660 mm)
- Load height: 31" (787 mm) above finished floor

Why purchase Auxo Medical Reconditioned Biomedical Equipment?

Auxo thoroughly inspects and re-certifies all of its reconditioned biomedical equipment. Batteries are checked or replaced, and each unit is cosmetically repaired. Calibrations are made and all safety inspections will be completed. In most cases our equipment warranties meets or exceeds those of the OEM.