



## Streamline 4 Pain Table - Fully Refurbished

Designed for high patient throughput for surgery centers, pain management, and other multidisciplinary medical facilities, the Streamline series X-ray imaging tables are the lowest priced tables in the imaging table industry. With a 24" top and carbon fiber diving board style layout, the streamline x-ray imaging table eliminates the need for arm boards and gives easy C-arm access into one innovative design. The operating table's light weight construction and 500 lb. weight capacity gives users easy transportation with some of their heaviest patients.

### Features:

- Carbon fiber, diving board style for easy C-Arm access
- 500 lb. patient capacity
- 84" x 24" table top
- 2" thick pad with facial cutout
- 48 Month Warranty
- 24" top eliminates need for arm boards
- Simple to use top quality wheel locks
- Easy patient transfer to table
- Hand control: activates all motions
- Patient straps: 2



### Specs:

- Height (28" - 38")
- Lateral Tilt (0° - 17°)
- Trendelenburg (0° - 17°)
- Longitudinal Travel (10° +/-)
- Length: 84"
- Width: 24"
- Weight: 590 lbs

### Why an Auxo Medical Refurbished Table or Stretcher?

Auxo tables and stretchers are refurbished to the some of the highest standards in the industry.

We completely disassemble each table and stretcher, and completely rework or replace all worn parts and return each table and stretcher to OEM standards. Each table and stretcher is reassembled, tested and is warranted against defects.